

RAIZER[®]

by liftup



Technical Specifications

- Lifting capacity: Max 150kg
- Lifting time: 20 - 30 seconds
 - Weight: Total kg
- Charger (2pieces) - Soneil 12V, 100-140V and 12V car adaptor
- Number of lifts when fully charged - 40 with max load and approx with average load.
 - Charging time: Max Hours

State of the Art Lifting Chair

Liftup's new and patented emergency lifting chair is a fast and safe solution when a person has fallen and needs help with transfer to either a standing or sitting position.

Benefits

- Mobile lifting chair for the transfer from floor to a sitting or standing position.
- Minimal physical effort for the assistant
- Both the fallen person and the assistant can feel comfortable and safe in the process
- Easy to transport in 2 light bags
- Assembly prior to operation: max 3 minutes
- Cleaning - friendly
- CE - Marked
- Patent Pending



HMS Medical

Specialising in patient transfers

T: 1300 106 003

F: 07 3036 6625

E: info@hmsmedical.com.au

www.hmsmedical.com.au