

## Which HoverMatt for What Patient

### Heat-Sealed

Heat-sealed construction eliminates needle holes that are potential bacteria pathways.

### Double Coated

The Double Coated version offers fluid and stain resistance, making it ideal for the OR or Labor & Delivery.

Heat Sealed and Double Coated are available in 4 size widths:

- (28") 71cm W x 198cm L
- (34") 86cm W x 198cm L
- (39") 99cm W x 198cm L
- (50") 127cm W x 198cm L

### Single Patient Use (SPU HoverMatts)

Works just like the reusable HoverMatt air transfer system but addresses infection control and reprocessing concerns.

SPU HoverMatts are available in 3 size widths:

- (34") 86cm W x 198cm L
- (39") 99cm W x 198cm L
- (50") 127cm W x 198cm L

### Which HoverMatt do I choose?

The weight and height of your patient / girth and the bed or trolley will determine the size mat you choose.

A standard bed accommodates a 34" matt and some will take a 39" matt. A hovermatt should not hang over the side of the mattress. as this reduces the effectiveness and in some cases can be dangerous.

A Bariatric bed is suitable for the 50" hovermatt.