

Bariatric Wedge

Consistently and reliably supports a comfortable side lying position for heavy weighted patients.

The unique design utilises an extra-firm core that supports heavier weights and a memory foam topper that optimises pressure redistribution, whilst the non-slip backing ensures it stays in position.

Features:

- Highly reliable turning solution for Bariatric / Plus size patients
- Manufactured of a very firm, high density core foam which will hold the torso in a turned position, while resisting the collapsing forces of heavier patient weights
- Layer of high-density Visco-Elastic Foam added to the patient contact surface, which enhances comfort and pressure reduction where it is needed most
- Ideal way to enhance your pressure management protocols; great for patients with pressure sores
- Non-slip fabric on the bottom of the reusable version prevents the wedge from migrating away from the patient
- This cushion maintains the patient's body position in bed, chair or other piece of equipment



Options Available:



BA1C-PFW001 (30cm)
30cm x 30cm x 20cm



BA1C-PFW002 (60cm)
60cm x 30cm x 20cm



BA1C-PFW003 (90cm)
90cm x 30cm x 20cm

PRODUCT CODE	PRODUCT NAME
BA1C-PFW001	Bariatric Wedge - 30cm
BA1C-PFW002	Bariatric Wedge - 60cm
BA1C-PFW003	Bariatric Wedge - 90cm

Distributed by:



12IILFH
7
)
(LPHVPHF)

12IILFH
7
)
(LPHVPHF)