

Patient Transfer Boards



Lateral Folding Transfer Board

- Is beneficial for the transfer of patients who find it difficult to lie flat or suffer from breathing difficulties.
 - 150cm x 74cm
 - Anti-slip pads ensure maximum grip during transfer
 - Six ergonomically placed handholds
 - High strength polypropylene for maximum durability
 - Easily wiped clean using soap & water, disinfectants or alcohol-based solutions (up to 75%)
- CODE: BAIC-QLTS**



These boards are designed to assist in the lateral transfer of patients and are NOT intended for lifting. They should only be used to bridge a gap of 15cm max between the surfaces and therefore they do not carry a maximum safe working load limit.

SallyBoard Transfer Boards (Slide Boards)

Lightweight for easy patient transfer.

- Constructed from lightweight, HDPE – much lighter than many transfer boards on the market
- Hanging hole at one end for tidy storage and easy accessibility
- 3 handles on each side for ease of manoeuvring
- Smooth top surface and curved edges for easy patient transfer
- 4 grip strips on underside to keep board in place
- Good level of flexibility in board where bed / barouche heights are not always at the same level
- Radiolucent – can be used for x-rays
- Large and medium size

Large :


180cm x 58cm W x 0.5cm Thick

CODE: MT-SBL

Medium:

150cm x 58cm W x 0.5cm Thick

CODE: MT-SBM



HMS Medical

Specialising in Patient Care

www.hmsmedical.com.au

T: 1300 106 003

Dir: 07 3117 9539

F: 07 3036 6625

E: sales@hmsmedical.com.au