

Bariatric Wedges



Bariatric Wedges provide a solution when turning and positioning Bariatric Patients. They are made using high density foam with a memory foam upper to give optimal comfort but firm to support the increased weight.


- The Convex shaped surface of the wedge makes it more supportive than standard turning and positioning wedges.
- Wedges have a non-slip base to prevent them moving when in use.
- The extra firm foam supports more weight and the memory foam top layer enhances comfort and pressure redistribution.
- Covered in a Polyurethane premium barrier textile that has superior stretch and breathability properties. Is skin friendly and shear reducing.
- Waterproof
- Covering is Anti-microbial and Anti-fungal
- Fire Resistant
- Latex free
- Seams are welded and sewn for added strength and to prevent possible ingress of any fluids.
- Available in 3 sizes



BARW9030



BARW6030



BARW3030

Code	Dimensions	Qty
BARW3030	30cm x 30cm x 20cm H	Each
BARW6030	60cm x 30cm x 20cm H	Each
BARW9030	90cm x 30cm x 20cm H	Each

HMS Medical

Specialising in Patient Care

T: 1300 106 003
Dir: 07 3117 9539
F: 07 3036 6625

E: sales@hmsmedical.com.au

www.hmsmedical.com.au