

Delta Body Bags

Light Weight

Material

Delta Medical body bags are a non-porous woven polyethylene fabric offering superior quality and is leak resistant, preventing leakage of body fluids while transporting and storing human bodies.

Features

- Body bags are not designed to be washed or re-used. They are discarded and incinerated after one use.
- The “D” shaped zipper system allows for easy access.
- Sewn ends
- Delta Medical body bags are manufactured in Brisbane and distributed nationally.

KEY SPECIFICATION

Product code	6003 – Baby 6002 – Child 6001 – Adult 6006 – Bariatric
Sizing	750x400mm – Baby 900x560mm – Child 2100x900mm – Adult 2400x2000mm – Bariatric
Quantity	1pcs/polybag, 20pcs/carton - Baby/Child 1pcs/polybag, 10pcs/carton - Adult 1pcs/polybag, 5pcs/carton - Bariatric
Material	White Poly Weave

----- Zipper opening

