

For easy identification of veins in the Newborn



The Neo Vista is a portable vein illumination device for neonates and babies.

When activated it displays a map of the vasculature on the surface of the skin, allowing clinicians to verify patency and avoid valves or bifurcations.

It is a more effective venepuncture procedure, with less patient discomfort, resulting in higher patient satisfaction and ensuring shorter procedure time.

- Easy to use, economical and hands free
- Increases the success rate with first attempt of venepuncture
- Splash proof, so that the user can clean it with disinfectant
- Minimum maintenance required

## SIZES AVAILABLE

When ordering please quote the following (minimum order qty 1 box):

ORDER#	DIMENSIONS	QTY/box
TRANS1	108mm x 32mm (max) x 38mm (max)	20

Distributed By:



T: 1300 106 003

07 3117 9539

F: 07 3036 6625

E: sales@hmsmedical.com.au

Also available in a Plus size.