

Air Supply Trolley

Hovertech HTAIR



AirSupply not included



- Stainless Steel Trolley that securely holds the Hovertech HTAIR AirSupply Motor
- Side sections to accommodate the Reusable or SPU HoverMatts and optional Coversheets
- Hooks to coil power cord around when not in use
- Black nylon base with 5 castors
- Made In Australia

Code: HMAST13



HMS Medical

Specialising in Patient Care

T: 1300 106 003
Dir: 07 3117 9539
F: 07 3036 6625

E: sales@hmsmedical.com.au

www.hmsmedical.com.au