

IndeeLift®

The People Picker Upper!

PPU-S For Home and Care Centers

Consumer Model

The "People Picker Upper" standing model (PPU-S) safely lifts a fallen person from the floor, eliminating the need for a call to emergency services while transforming a stressful fall into a non-event.

It can raise individuals up to 181kg to a seated or standing position, and allows for easy transfers to and from a wheelchair, bed, commode, exercise equipment or even the floor.

The standing lift function allows users to be safely raised to access a high bed or an overhead cabinet - a step stool without the step!

Can be self-operated or with the assistance of a caregiver.

Restores dignity and independence.




Assisted Operation



Self Operation




**Standing
Function**



Summary of Features

- Beveled seat-front for easy slide on transfer surface
- Wired remote for convenient operation
- Swinging/Removable Rise-Assist Handles for easy secure transfer
- Dual safety belts for added safety
- Powder coated for easy cleaning
- Small footprint for tight spaces
- Lithium-Ion battery powered with automatic charger
- 9' (2.74m) power cord



Physical Specifications

- Safe Working Load 181 Kg
- Overall Depth 26.5" (67 cm)
- Maximum Overall Height (raised) 70.75" (174 cm)
- Minimum Overall Height (lowered) ... 41.75" (100 cm)
- External Width incl Rise Handles..... 23" (58 cm)
- Maximum Seat Height 30" (76 cm)
- Seat Width 14" (35 cm)
- Wheel Diameter..... 5" (13 cm)
- Unit Weight 78 lbs. (35 kg)

Electrical Specifications

- Standard 2-prong 110 VAC Power Cord (USA)
- Standard 2-prong 220 VAC Power Cord (International)
- Operating Environment: >+5° C <+40° C (41° to 104°F)
- Lithium-Ion Battery: 24 VDC

Maintains Privacy, Dignity and Independence



T 1300 106 003
Dir: 07 3117 9539

E: sales@hmsmedical.com.au

www.hmsmedical.com.au