

# Tube Trees

Anaesthetic Circuit and Line organiser


HMS Medical

Specialising in patient transfers



180 degree lateral adjustment

Vertical Adjustment



Code: HMS-04



Code: HMS-03

- Easily secures anaesthetic tubes, suction, gas sampling lines
- Manages adult and paediatric corrugated circuits
- Slots are tapered to hold the circuit tubing securely
- Made from 4.5mm acrylic which can be easily wiped down using a neutral PH solution
- Safe for MRI use
- Available in static or articulating styles
- Base slides under the trolley or bed mattress
- Height 31cm x Width 12.5cm x Depth 21.5cm

Dir: 07 3117 9539

T: 1300 106 003

F: 07 3036 6625

E: sales@hmsmedical.com.au

Latex Free  
Made in Australia

[www.hmsmedical.com.au](http://www.hmsmedical.com.au)

HMS-TTBroch-032019