

IGuard


Patient Eye Protectors



LATEX FREE | SINGLE USE | DISPOSABLE



- Rigid, clear protective patient eye covers.
- Designed for use in Operating Rooms, Emergency Departments, and during general procedures where there is a risk of corneal or orbital damage from foreign bodies, splash or pressure.
- Non-allergenic self-adhesive foam cushion for easy, accurate application and removal
- Stays in place when patient is supine, prone or lateral
- Latex free
- Non Sterile
- Single Use
- No micropore or tape required



Code	Description	Qty
SM 9-0210-001/bag	IGuards, Adult eye protectors. Packs of 25	Bag of 25
SM 9-0210-001/box	IGuards, Adult eye protectors. 10 packs per box	Box of 250

T: 07 3117 9539
1300 106 003
F: 07 3036 6625

E: info@hmsmedical.com.au

www.hmsmedical.com.au