

Ergonomic Seating

Chairs & Stools

ERGO STOOL - with Black Base

- Fully upholstered Seat in in high grade Black Vinyl.
- Moulded seat foam for maximum comfort and durability.
- 480mm Span Black base is both compact and stable.
- 160 Stroke Gas Lift allows for excellent range of height adjustment.
- Circular ring Gas activation allows for easy access via either hand.
 - Seat Height: 520mm - 720mm
 - Seat Dimension: 390mm x 90mm



CODE: HMS ES

ERGO STOOL - with Aluminium Base and Strut

- Fully upholstered Seat in high grade Black Vinyl.
- Molded seat foam for maximum comfort and durability.
- 500mm Span Polished Aluminium base is both compact and stable.
- 200mm Stroke Gas Lift allows for excellent range of height adjustment.
- Circular ring Gas activation allows for easy access via either hand.
 - Seat Height: 520mm - 720mm
 - Seat Dimension: 390mm x 90mm



CODE: HMS-ES-AL

PU STOOL

- Simple, functional design - easy to operate.
- Seat made from high quality polyurethane.
 - Hard wearing and easy to clean.
- 480mm Span Black base is both compact and stable.
- 160 Stroke Gas Lift allows for excellent range of height adjustment.
- Circular ring Gas activation allows for easy access with both hands.
 - Seat Height: 440mm - 580mm



CODE: PU Stool

ERGO DRAFTING STOOL with Footring

- White upholstered seat in high grade vinyl.
- Moulded seat foam for maximum comfort and durability.
 - Kit comes with 16" adjustable footring.
 - 500mm span polished aluminium base.
- 270mm stroke gas lift allows for an excellent range of height adjustment.
- Circular ring gas activation allows for easy access with either hand.
 - Seat Height: 580mm - 830mm



CODE: HMS ESTall

T: 1300 106 003

Dir: 07 3117 9539

F: 07 3036 6625

E: sales@hmsmedical.com.au

HMS Medical

Specialising in Patient Care

www.hmsmedical.com.au

Ergonomic Seating

Chairs & Stools

ERGO STOOL WITH BACK SUPPORT

- Upholstered in high grade white or black vinyl.
- Firm cushion seat and supportive wing back.
- Matched with 2 lever semi ergo mechanism.
- White has polished base and chrome gas lift and clear PU castors.
- Black has a black nylon base with black gas lift and twin wheel castors.

Code: HMS BES-W/B (Black Seat)

HMSWES-W/B (White Seat)



PU STANDARD CHAIR

- Seat and Back made from High quality polyurethane.
 - Hard wearing and easy to clean.
 - 2 Lever ergo mechanisms.
- Black nylon base with black gas lift and twin wheel castors
 - Seat Height 415mm - 545mm

Code: HMS-PU-STD



PUTALL CHAIR

- Seat and Back made from High quality polyurethane.
 - Hard wearing and easy to clean.
 - 3 Lever ergo mechanism.
- Black nylon base with black gas lift and twin wheel castors.
 - Drafting height gas lift and Footring.
 - Seat Height 580mm - 840mm.

Code: HMS-PU-Tall



HMS Medical

Specialising in Patient Care

T: 1300 106 003

Dir: 07 3117 9539

F: 07 3036 6625

E: sales@hmsmedical.com.au

www.hmsmedical.com.au

Ergonomic Seating

Chairs & Stools

SADDLE CHAIR WITH BACKREST

- Creates correct seating posture to assist with back complaints and blood flow.
 - Seat and back upholstered in high grade vinyl.
 - Adjustable backrest.
 - Matched with 3 lever fully ergo mechanism.
- Black nylon base with black gas lift and twin wheel castors.
 - Seat Height 530mm - 650mm

Code: HMS Saddle W/Back



SADDLE STOOL

- Creates correct seating posture to assist with back complaints and blood flow.
 - Seat upholstered in high grade vinyl.
- Black nylon base with black gas lift and twin wheel castors.
 - Seat Height 530mm - 650mm

Code: HMS Saddle



PU SIT STAND STOOL

- Seat and back made from high quality polyurethane.
 - Tilt and height adjustable.
- Black nylon base with black gas lift and twin wheel castors.
 - Seat Height 580mm - 840mm

Code: HMS PU SSS



PU HAND HOLD STOOL

- Seat made from high quality polyurethane.
 - Convenient hand hold grip.
- 160 mm Stroke Gas Lift allows for an excellent range of height adjustment.
- Round gas lift activation ring allows height adjustment by either hand.
 - Seat Height 440mm - 590mm

Code: HMS PU HHS



HMS MedicalSM

Specialising in Patient Care

T: 1300 106 003

Dir: 07 3117 9539

F: 07 3036 6625

E: sales@hmsmedical.com.au

www.hmsmedical.com.au