

Replacement Pads

Replacement Mattresses, Armboard Pads and YelloFin stirrup inserts



- Made in Australia
- Multi layer foam for comfort and protection of the patient.
 - Seams sealed to prevent fluid ingress
 - 2 way stretch fabric
- Mattresses with or without perineal cutout
- Different Mattress configurations available

Code	Dimensions
------	------------

Mattress


ST-I103-60	3 Piece 60mm thickness
ST-I104-60	4 Piece 60mm thickness
ST-I103-80	3 Piece 80mm thickness
ST-I104-80	4 Piece 80mm thickness
ST-I103-100	3 Piece 100mm thickness
ST-I104-100	4 Piece 100mm thickness

Armboard Pads

ST-I03	200mm wide x 60mm thick
ST-I03-I660	160mm wide x 60mm thick
ST-I03-I680	160mm wide x 80mm thick
ST-I03-I610	160mm wide x 100mm thick

YelloFin Stirrup Inserts

YF01	Replacement inserts for YelloFins
------	-----------------------------------



Code	Dimensions
------	------------

Positioning Wedges

ST-300	300mm W x 130mm H x 300mm L
ST-375	375mm W x 170mm H x 375mm L
ST-560	560mm W x 425mm H x 380mm L
ST-M300	300mm W x 130mm H x 380mm L