

Safety Straps for Patient Positioning

SPU Patient Positioning Straps



Disposable patient positioning straps are used to assist in immobilizing patients and are great for infection control purposes. Proper patient positioning is critical to patient safety, a durable safety strap is often required for specific positioning techniques and positions. Our SPU Positioning Straps will reduce the risk of patient injury and maximise patient safety. They are made from a soft flexible material with a strong velcro closure.

- The velcro closure gives versatility to size
- Not rail specific
- No metal clips or closures.
- SPU - no infection control concerns
- Latex free
- 2100mm x 160mm wide.
- Boxes of 36
- Boxes of 6 available

SBW2100 - Box of 36
SBW2100 -6 Box of 6



Reusable Patient Positioning

A Reusable range of positioning straps are also available.

Velcro closure still gives the versatility of size.

Are made from a soft Carflex Material which can be wiped over with neutral PH cleaning solution.

- Comfortable and reduces risk of patient injury.
- Available in 2100mm / 2500mm / 3000mm x 150mm wide
- Not rail specific
- No metal clips or closures.
- Latex free

STR210
STR250
STR300



Direct: 07 3117 9539

T: 1300 106 003

F: 07 3036 6625

E: sales@hmsmedical.com.au

HMS Medical

Specialising in Patient Care

www.hmsmedical.com.au