



# SYST'AM® P311CD / DUOGEL CUSHION

GEL CUSHION MADE OF VISCOELASTIC GEL AND MULTIBEARING SURFACE

## MATERIALS



Visco Gel



Removable Polymaille® cover



NF EN ISO 597 - 1 & 2

Cleaning of the gel :



(Do not immerse in water)

Cleaning of the cover :



(UP TO 90°C)

## INDICATIONS

### Risk



P311CD / DUOGEL-Cushion can be used, when :

- Seated position is necessary resulting from diseases or handicap and there is at the same time an increased risk of pressure sore.

SYST'AM® P311CD  
(without cover)



SYST'AM® P311CD  
(with cover)





| Item code     | Designation              | Sizes (cm / inches)<br>(W x L) | Patient weight (kg / lbs) |           |
|---------------|--------------------------|--------------------------------|---------------------------|-----------|
|               |                          |                                | Min.                      | Max.      |
| P311CD42421HW | P311CD - DUOGEL® cushion | 42 x 42 / 17 x 17              | 40 / 88                   | 140 / 308 |

## PRODUCT «ADVANTAGES»

### VISCOELASTIC HIGH-RETICULATED POLYURETHANE GEL

- Avoids gel migration and prevents the bottoming effect of the bony protuberances.



### 100% IMPERMEABLE

- Completely adapted to incontinent patients.

### AN INNOVATIVE MULTI-BEARING SURFACE

- With an innovative technology, the Duogel reduces greatly the maximum pressure exerted on sensitive areas by increasing the contact surface between the patient and the cushion.

#### WHITE AREA :

A low bearing, viscoelastic, self-contained and leak-free polyurethane gel for minimum pressure on high-risk areas such as pelvic protuberances.

#### BLUE AREA :

A high bearing surface reticulated viscoelastic P.U. gel on low risk areas for improved comfort and patient stability.

Distributed By:



T: 1300 106 003  
Dir: 07 3117 9539  
F: 07 3036 6625  
E: sales@hmsmedical.com.au

