



The Neo-Guard Family has been designed to protect the sensitive nasal and facial skin of the neonate undergoing oxygen treatment. They guard the neonatal septum from trauma and irritation which is often caused by pressure & motion from the oxygen prongs and masks.

- Hydrocolloid material is very skin friendly even to the sensitive skin of neonates.
- Repositioning of the Neo-Guard on neonates is easy and does not cause any discomfort.
- The special hook & loop design at the bottom of every piece of Neo-Guard provides an extra method of securing multiple types and sizes of nasal device.
- Easy observation of the condition of the nasal nares.
- Multiple sizes & designs available for the best match of therapies to baby weight.

Optimization

Three unique designs to best suit the type of oxygenation therapy.

Neo-Guard

- Specially designed for nasal CPAP
- It adopts a partially cut-through star around the holes to generate retention on the prongs, effectively securing them in position and providing a better seal for stable alveola pressure.

Neo-Guard Plus

- Intergrates a layer of foam with a silicone base.
- Will give better protection against potential pressure necrosis to the columella.



Code	Dimensions
NG000	CPAP Neo-Guard <700g
NG001	CPAP Neo-Guard 700g - 1250g
NG002	CPAP Neo-Guard 1250g - 2000g
NG003	CPAP Neo-Guard 2000g - 3000g
NG004	CPAP Neo-Guard >3000g
NGP-S-M	CPAP Neo-Guard w/foam <2000g
NGP-M-L	CPAP Neo-Guard w/foam >2000g
HFNG-S	High flow Neo-Guard <1250g
HFNG-M	High Flow Neo-Guard 1250 - 2000g
HFNG-L	High Flow Neo-Guard >2000g

High Flow Neo-Guard

- The holes are significantly larger than the cannula, providing a clear flow path.
- Contact between the nasal nares and the cannula is minimal therefore reducing any potential trauma or irritation.

