



# Neo-Shield

- Latest technology against heat radiation and conduction

- silver reflective surface
- insulation foam layer

Providing accurate temperature monitoring.

- Silicone based adhesive: the most gentle adhesive in modern medicine, especially for neonatal skin.
- Maintains skin integrity; does not cause pain or damage when removed.
- Will not leave residue on the skin.
- Able to be repositioned without losing adhesiveness:  
Extends product life; cost effective.
- Now available in 2 sizes

When ordering please quote the following

Description	Order #	QTY
Neo-Shield 2.5cm x 2.5cm	REF-001	50 pieces / Box
Neo-Shield 3.5cm x 3.5cm	REF-002	50 pieces / Box



## HMS Medical

*Specialising in patient transfers*

T: 1300 106 003

Direct: 07 3117 9539

F: 07 3036 6625

E: [sales@hmsmedical.com.au](mailto:sales@hmsmedical.com.au)

[www.hmsmedical.com.au](http://www.hmsmedical.com.au)

