E: sales@hmsmedical.com.au

M2 300 kg Wide El-Tip

The shower commode chair for the largest users



Big advantages with M2 300 kg Wide El-Tip



- M2 300 kg Wide El-Tip can carry users with a weight up to 300 kg.
- With 72 cm between the armrests and an extra wide backrest, most users can be placed in this chair.
- The chair is equipped with wheels with double action brakes. The diameter is increased from 125 mm to 150 mm which makes it very easy to maneuver.
- It has a weight of 32 kg and yet a simple construction that is easy to clean.
- Reinforced armrest and soft foam seat.
- With electric tilt function, the chair can be tilted 8° forward so that users who is capable of standing can be helped into the chair and out again.
- For lift users it can be tilted up to 28° backwards, so it's easy to get the user into the chair. That way the staff's working environment is ensured as well.
- The seat height is 61 cm, while the negative tip only has an entrance height of 53 cm.
- As with all our products, it is possible to customize the chair according to the needs of the user either through accessories or special alterations.







Technical data

| M2 300 kg Wide El-Tip | | |
|---------------------------------|---|----------|
| Seat height | | 61 cm |
| Seat width | 75 cm | |
| Seat depth | 46 cm | |
| Max. width | 83 cm | |
| Space between armrests | 72 cm | |
| Length (without the footrests) | 83 cm | |
| Weight | 32 kg | |
| Max. load | 300 kg | |
| Seat angle in degrees | 28° backwards, 8° forward | |
| Entry height | | 53 cm |
| Wheels | 150 mm wheels w/grey rubber and ball bearings. | |
| | Double-acting brakes on 3 wheels. | |
| | 1 directional fixed wheel. | |
| Height adjustable | No 20 v 70 am | |
| Backrest | 30 x 70 cm. | |
| | Depth adjustment up to 10 cm. | |
| | The net fabric with sewn did and 4 | |
| | buckles for depth adjustment. | |
| Foam parts | Seat: 5,5 cm foam | |
| | Armrests: 2,5 cm PU foam | |
| Electric motor / Battery system | Stepless blocked actuator. | |
| | Batteries can usually drive the chair 50+ times back and forth at | |
| | full charge. EL-tightness = IP66 (Dust tight and protected against water beams) | |
| Surface / metal | Stainless steel. Powdercoated. | |
| Product | Art. no. | ISO. Nr. |
| M2 300 kg Wide El-Tip | 310218 | 09 12 03 |



T: 1300 106 003

07 3117 9539 F: 07 3036 6625

E: sales@hmsmedical.com.au



Midtager 22 - 2605 - Brøndby - Denmark

Tlf.: +45 4396 2400 - E-mail: info@hmn.dk - Web: www.hmn.dk

HMN a/s

