

Zero pressure on the Occipital area helps to prevent decubitus ulcers, pressure alopecia and fluid pooling

Fingered lateral support provide excellent pressure distribution and stable head support with comfortable fit

Raised neck area supporting cervical vertebrae and trapezoid muscle

Smooth contoured surfaces allows maximum patient comfort

Non-slip base grips to operating table mattress

Concave base aids self adjusting motion and allows air ventilation

Air ventilation through the base maintains normothermia at the occipita bone

VSHHP23 Supine Head Positioner Medium Style 2
VSHHP24 Supine Head Positioner Large Style 2

Unique and innovative design that **redistributes weight of the head** and provides enhanced pressure relief compared to regular Head Rings.

The Positioners are made from two components. The inside of the product is comfortable Visco-elastic foam. The outer coating is a seamless, non-porous **easy to clean Hi-Tech Membrane with high tear strength** making it a very robust product.

Suitable for **Robotic, ENT, Ophthalmic and General Surgery** in the Supine Position.

VSKP13 Supine Knee Supine Medium (Pair)

Unique and innovative design with a proven track record that **supports the Leg and Knee** and provides enhanced pressure relief under the Knee compared to regular Patient Positioners.



Zero pressure under the knee area helps to prevent pressure sores, pressure build up around the popliteal artery and peripheral ligaments

Raised knee area supporting leg and preventing hyper-extension of the knee joint

Raised exterior sides help to prevent leg rolling or movement from the leg when in a relaxed state

Concave base aids self adjusting motion and allows air ventilation

Non-slip base grips to operating table mattress

VSFP13 Supine Foot Positioner Medium (Pair)

Unique and innovative design that **redistributes weight of the Foot** and provides enhanced pressure relief under the **Calcaneus Bone** compared to regular Heel Pads. strength making it a very robust product.

The **easy to clean** outer coating has **excellent Puncture/Burst Resistance**. It has excellent flexibility that **improves conformity to the Foot**.

The Supine Knee Positioners and Foot Positioners are intended to be used simultaneously.

Two Year Product Warranty

Zero pressure on the calcaneus bone helps prevent decubitus ulcers developing

Raised sides protects the foot from lateral force and leg rotation

Smooth contoured surface allows maximum surface redistribution

Unique non-slip base grips to operating table mattress

Ergonomical support surface contoured to the back of the leg and the achilles tendon

Concave base aids self adjustment motion and allows air ventilation

Air ventilation through the base allows normothermia at the calcaneus bone



VascoCare Medical Limited
Enterprise Research Incubation Centre
South East Technological University
Kilkenny Road, Carlow
Ireland
T: +353 (0)86 3513127
e: customercare@vascocare.com
w: www.vasvocare.com

Unique Code: PI-VSP-01 Date of Issue: 9th Nov 2022 Revision: 0
© The information provided is the copyright of VascoCare Medical Limited.

EVOLUTION PATIENT POSITIONERS

EVOLUTION Supine Positioners

| | | |
|------------|---------------|--|
| REF | VSHP23 | Supine Head Positioner Medium Style 2 |
| REF | VSHP24 | Supine Head Positioner Large Style 2 |
| REF | VSKP13 | Supine Knee Positioner Medium (Pair) |
| REF | VSFP13 | Supine Foot Positioner Medium (Pair) |

Intended Use

These Reusable Patient Positioners are used to protect the bony prominence of a Patient's body and help reduce the risk of developing pressure sores while a Patient is resting on the Operating Table or other Support Surface during Medical/Surgical Procedures and Examinations by providing an ergonomically shaped stable support to hold the Patient and preventing them from accidentally moving.



Supine Foot Positioner



Supine Head Positioner

Supine Knee Positioner

VascoCare Medical introduced innovative Supine Positioners to the EVOLUTION Patient Positioner Range. These provide enhanced pressure relief with a proven anatomical design.

| Product Description | Length (±5mm) | Width (±5mm) | Height (±5mm) | Weight (±20g) |
|---------------------------------------|---------------|--------------|---------------|---------------|
| Supine Head Positioner Medium Style 2 | 250 | 210 | 116 | 235 |
| Supine Head Positioner Large Style 2 | 282 | 237 | 132 | 360 |
| Supine Knee Positioner Medium (Pair) | 298 | 382* | 141 | 560 |
| Supine Foot Positioner Medium (Pair) | 261 | 366* | 130 | 700 |

* Dimension when pair of Positioners are side by side.

Features

- Unique and innovative design that **redistributes weight of the patient's body and provides enhanced pressure relief** compared to regular Patient Positioners.
- The positioners are made from two components. The inside of the product is comfortable **Visco-elastic foam**. The outer coating is a **seamless, non-porous Dual Coated Hi-Tech Membrane** with high tear strength making it a very robust product.
- The outer coating has **Excellent Abrasion Resistance** and is **Puncture/Burst Resistant**.
- It has excellent flexibility that **improves conformity to the skin**.
- Provides **Zero Interface Pressure** under the Bony Prominences.
- High strength bond between foam and coating creates an ultra strong laminated structure for a **long life product**.
- The outer coating provides an excellent **fluid-proof barrier without any seams or other bacteria harbouring features on the outer surface**.
- Completely **Latex free**.
- **Compatible with CT Scanners** as these do not cause artefacts.
- Head Positioners provide **excellent Neck Support and Head Stability** compared to regular Head Rings.
- Knee Positioners provide **excellent Support under the Leg** compared to regular Domed Patient Positioners.
- Foot Positioners provide **excellent Foot Support with Zero Pressure under the Heel** compared to regular Heel Pads.
- Competitively priced to make it a **cost effective option** against regular positioners.
- **Two Year Product Warranty**.

Distributed By



T: 1300 106 003

Dir: 07 3117 9539

E: sales@hmsmedical.com.au

