

Body Positioners for Colonoscopy

Lateral positioner to keep the patient safely on their side.

- Bevelled edges to make sliding under the mattress easier.
- 8mm thick acrylic.
- 500mm wide x 500mm high x 400mm deep.
- Comfortable foam pad with welded seams.
- Latex free
- Available with or without pad.
- Replacement pads available.
- Easy to clean - wipe over with standard cleaning wipes.
- Do not submerge in water, autoclave or place in washer.



Reusable Patient Positioning straps

- Velcro closure - gives versatility of size.
- Made from soft Carflex material for comfort.
- Not rail specific - easily applied and easily adjusted.
- Available in 3 lengths - all 150mm wide.
- No metal clips or buckles
- Latex free
- Wipe over with standard hospital cleaning wipes.
- Do not submerge in water, autoclave or put in washer.



Code	Description / Size	Qty
HMS05	Complete with pad	Each
HMS06	Lateral support frame only	Each
HMS07	Pad only for HMS06	Each
STR300	Reusable positioning belt. 3000mm x 150mm	Each
STR250	Reusable positioning belt. 2500mm x 150mm	Each
STR210	Reusable positioning belt. 2100mm x 150mm	Each

HMS Medical

Specialising in Patient Care

T: 1300 106 003
Dir: 07 3117 9539
F: 07 3036 6625

E: sales@hmsmedical.com.au

www.hmsmedical.com.au