



Isolux Headlight / proCam combination improves operative site visibility and video capture



proCam 4K MicroCamera



proCam 4K MicroCamera

World-leading camera to improve your everyday communication and documentation

proCam 4k Headlight Camera - first miniature 4k camera. Features a 13 megapixel sensor that captures both 4k video and high resolution still photos from proCam / headlight combination. The proCam's powerful sensor allows deep, clear magnification. It's easy to capture stable, super high quality video and images any time with no interruptions to the procedure.

All advantages in one camera. Futudent non-interruptive cameras cover most of the demands of documentation: full HD videos and high quality stills (up to 13 mb). Cloud services make it fast and secure to share videos and images to colleagues and clients.

Key Features

- 4k Video Quality
- Live Streaming, Recording and Still Photo
- Headlight or Fixed Mount
- Fixed Focus

Applications

- Clinical Education
- Live Surgical Demonstration
- Virtual Surgical Support

Camera warranty: 2 years standard

CAMERA SPECIFICATIONS

Resolution / Video	4K UHD, 30 fps (video)
Resolution / Still	8 / 13 mpx
Size (w/o lens)	32 x 28 x 23 mm
Weight (body only)	23 grams
Built-in zoom	3x
Strength	UHD video, high resolution, documentation
Connection	USB 2
Lens Options	<ul style="list-style-type: none"> • Standard 16mm • Hi magnification 25mm UHD • Wide angle 12mm and 8mm
Compatible with	Windows, Mac & Android

futudent»
SHARE your CARE

Futudent non-interruptive cameras are designed for taking great still photos and videos. There are dozens of applications for video in today's clinical setting and medical education.