

Arm Board Reusable Positioning Straps

The Reusable Arm Board Strap holds the patient's arm securely on the armboard during surgery and will fit most armboards. Made from Carrflex Extreme which is a soft, flexible and breathable fabric with velcro closure. Carrflex is Anti-microbial, Anti-fungal and fire resistant. Latex free.

900mm long x 100mm wide.

Code	Dimensions	Qty Buy
PAS9010S1603	900mm x 100mm	Eaches



Code	Dimensions
2103	25cm x 5cm x 1cm
2103-1	32cm x 5cm x 1cm



Arm / Hand Reusable Positioning Strap

The Gel positioning strap will keep the patients arm secure when out on an armboard. The gel pad is comfortable and will not cause damage to patients skin while the length is easily adjusted with the hook and loop closure.

HMS Medical

Specialising in Patient Care

T: 1300 106 003
Dir: 07 3117 9539
F: 07 3036 6625

E: sales@hmsmedical.com.au

www.hmsmedical.com.au